

In tasks with high motion collinearity, noisy ICs must be added to the regression model after orthogonalization to the signals of interests and other BOLD-related ICs.




Fig. 1. FD vs DVARS of the denoising pipelines for a representative subject. Each dot is a timepoint, each line a session. Optimal Combination of ME data improves motion denoising compared to single echo, ME-ICA is even more effective.




Figure 2. ICC(2,1) of the CVR maps from each pipeline.

In both figures:

- pre: before denoising
- echo-2: single echo data
- optcom: optimally combined (OC) ME data
- meica-agg: OC data, denoising with "noise" ICs timeseries
- meica-ort: OC data, denoising with "noise" ICs timeseries orthogonalised to CO2 trace
- meica-con: OC data, denoising with "noise" ICs timeseries orthogonalised to CO2 trace and "signal" ICs

ICA-based denoising strategies in highly motion correlated tasks with Multi-Echo BOLD fMRI

Stefano Moia¹, Maite Termenon¹, Eneko Uruñuela¹, Rachael Stickland², Molly Bright², César Caballero-Gaudes¹

1. Basque Center on Cognition, Brain and Language, Donostia, Spain, 2. Northwestern University, Chicago, IL



Northwestern

Background

- Functional MRI presents different sources of noise.
- Removal of task-correlated noise, such as in a Breath-Holding (BH) paradigms, is particularly challenging due to the tradeoff between effective denoising and signal conservation [1].
- Independent component analysis (ICA) can also be employed to identify noise and use this information during nuisance regression [2].
- Multi-echo fMRI allows increasing the contrast to noise ratio of the signal by optimal combination of the echoes, and identifying non-BOLD sources from BOLD-related components via ME-ICA [3-5].
- Main aim:** Evaluate different denoising variants to clean ME-fMRI data acquired during a BH task, in order to improve Cerebrovascular Reactivity (CVR) estimation.

Results

- Fig. 1: ME-fMRI with optcom and ME-ICA enhances denoising of motion-related effects (i.e. reduced FD-DVARS correlation), but the type of model with data-driven ICA nuisance regressors is key for denoising performance.
- Fig. 2: Optcom and meica-con yield the most reliable CVR maps, while meica-aggr is the least reliable since the effect of interest (CVR) is removed due to high collinearity of motion-related IC components.

Overall, meica-con produces the most reliable CVR maps with improved motion denoising.

More results in the companion website.

Bibliography

- Bright, M.G., et al. (2009), Neuroimage;
- Caballero-Gaudes, C. & Reynolds, R. (2017), Neuroimage;
- DuPre, E., et al. (2019), doi:-10.5281/ZENODO.2558498;
- Kundu, P., et al. (2012), Neuroimage;
- Kundu, P., et al. (2013), PNAS;
- Bright, M.G. & Murphy, K. (2013), Neuroimage;
- Moia, S., Stickland, R., et al. (2020), accepted in Proceedings of the 42nd Annual International Conference of the IEEE Engineering in Medicine and Biology Society (EMBC), Montreal, Canada;
- Power, J.D., et al. (2012), Neuroimage;
- Smyser, C.D., et al. (2011), Neuroimage;
- Chen, G., et al. (2018), HBM

Methods

- Seven healthy volunteers underwent 10 MRI sessions in a 3T Siemens PrismaFit scanner, spaced 1-week apart at the same time of day.
- A BH task adapted from [6] was administered at each session while collecting ME-fMRI data. CO2 levels were measured using a nasal cannula with gas analyzer (ADInstruments) and BIOPAC MP150 system. A T1-w image was collected during each session. The parameters can be found in the website version.
- ME data were decomposed using ICA (tedana [3]), and the ICs were manually classified into "signal" and "noise" ICs (see website).
- In order to obtain CVR maps, data preprocessing and analysis followed the steps described in [7] (optcom). The same steps were applied to the second echo volume (echo-2). Additionally, the optimally combined data were analysed using a similar pipeline, but including in the GLM the nuisance regressors: (I) the "noise" ICs timeseries (meica-aggr), or (II) the "noise" ICs timeseries orthogonalised w.r.t. the CO2 trace (meica-ort), or (III) the "noise" ICs timeseries orthogonalised w.r.t. the CO2 trace and the "signal" ICs timeseries (meica-con).
- FD and DVARS [8,9] were computed before realignment (pre) and on the optcom volume, after removal of the reconstructed volume from the nuisance regressors.
- ICC(2,1) was computed using 3dICC (AFNI) [10] on the CVR maps, obtained from the fit of the CO2 regressor in each of the data analysis pipelines.

Acknowledgements

Research supported by the Eunice Kennedy Shriver National Institute of Child Health and Human Development of the National Institutes of Health under award number K12HD073945, the European Union's Horizon 2020 research and innovation program (Marie Skłodowska-Curie grant agreement No. 713673), a fellowship from La Caixa Foundation (ID 100010434, fellowship code LCF/BQ/IN17/11620063), the Spanish Ministry of Economy and Competitiveness (Ramon y Cajal Fellowship, RYC-2017-21845), the Spanish State Research Agency (BCBL "Severo Ochoa" excellence accreditation, SEV-2015-490), the Basque Government (BERC 2018-2021 and PIBA_2019_104), the Spanish Ministry of Science, Innovation and Universities (MICINN; FJCI-2017-31814)